



Home HealthCare Laboratory of America
320 Premier Court
Suite 220
Franklin, TN. 37067
Service @hhla.com

LABORATORY UPDATE

HHLA's commitment to improve the quality of our service has led to a change in methodology and reagents for Drug testing.

The BIG change is the stability of the sample.

Refer to the following list:

- ★ Digoxin is stable for 3 Days
- ★ The following drugs are stable for 48 hours.

**Amikacin
Carbamazepine
Gentamicin
Phenobarbital
Phenytoin (Dilantin)
Theophylline
Tobramycin
Valproic Acid
Vancomycin**

Reference Ranges will be listed on the lab profile report with results. (No change)

Note: This change will mean HHLA will **NOT** be able to add on drug testing to samples after the stability time line.

If you need additional information or have any questions, please call Customer Service at 888-522-4452!

Effective immediately 06-06-2006